



SPECIALIST QUALITY CONTROLS

Exclusive worldwide manufacturers of Serum Indices Quality Controls

How confident are you that your laboratory results are accurate?

- Why use Serum Indices Quality Controls?
- Haemoglobin, Icterus, Lipaemia – 3 major interferences in laboratory analysis
- Millions of samples analysed every day at risk of being inaccurate
- Time wasted in re-sampling
- Time wasted in re-analysis
- Automated analysers automatically report Serum Indices results with NO quality control check

Use Serum Indices Quality Controls to remove risk of inaccuracies

Product features

- Liquid frozen product
- Haemoglobin, Bilirubin, Lipid, combined in one multi-parameter control
- 7 day open vial stability @ 2-8°C
- Options available for all Clinical Chemistry Analysers

Where can I buy?

Contact: C: **ADP Diagnostics**
Carrickfergus. Northern Ireland. U.K.

T: **+44 (0) 2893 366929**

E: **adpni@adp.ie**

W: **www.adpdiagnostics.com**

Ask us for details of our network of international distributors