


CareStart™ G6PD




Applications

- The CareStart™ G6PD Test is a visual screening test to identify G6PD (Glucose-6-phosphate dehydrogenase) deficient patients using whole blood sample.
- The CareStart™ G6PD Test can be a semi-quantitative assay in conjunction with a G6PD reader.

Features

- Rapid Diagnostic Test (RDT)
- Point-Of-Care Test with no specimen preparation
- Visual dye colorization method
- Fast and reliable test ($\geq 95\%$ sensitivity for samples defined by WHO as Class I and II G6PD deficiency)
- Simple and easy to use
- Results in 10 min.
- Broad range of assay temperature (18–32 °C)
- Store at room temperature
- CE Marked

Cat. No	Specimen	Specimen volume	Test time	Pack size
G0221	Whole blood	2 μ l	10 min	25 or 50 tests/ box



65 Clyde Rd., Suite A,
Somerset, NJ 08873
T. 732-873-4040
F. 732-873-4043
E-mail. info@accessbio.net
www.accessbio.net