

# SIMPLE OPERATION

Accuracy and simplicity



1. Collect Sample



2. Insert into Cuvette



3. Insert into Channel

## 3 Key Technologies INNOVATION

### Build-in Battery

- Charge 3Hr ,run 30Hr non-stop

### 3-level Calibration System

- Assure reliable and accurate results


### Multi-methodologies

- Reagent platform sharing

# DETAILS

Convenient,Fast, and Meticulous

## 3-level Calibration System



### TECHNICAL DATA

Dimensions		Applicable Methodologies
W × D × H	410*355*280mm	Biochemistry
Weight	5kg	Nephelometry
Operation		Turbidimetry
Power supply	AC 110~265 V	Immunofluorescence
Sample test time	< 5 minutes	Calibration
Patient Input	Bar code,RFID	Hardware calibration
Quality Control		Network calibration
Storage the newest reagents'curves updated by network connection to host server		Third party calibration

# TEST MENU

Multiple available reagents

Infection
PCT
Full range CRP
SAA
Renal-Function
HbA1C
CYS-C
mAlb
NGAL
RBP
β2-MG
α1-MG
Microelement
★ 25-OH-VD
Liver Function
★ ADA
TRF
PA
★ HPT

Cardiovascular disease
H-FABP
MYO
cTnI
HCY
NT-ProBNP
Lp(a)
Lp-PLA2
CK-MB
D-Dimer
Gastric Function
★ HP
PGI
PGII
Tumor Marker
★ ADA
AFP
CEA
★ PSA
★ fPSA

Reheumatism
ASO
CRP
Anti-CCP
RF
Autoimmunity
IgA
IgG
IgM
Hormones
HCG
★ LH
Thyroid Function
★ T3
★ T4
★ TSH
★ FT3
★ FT4
more to come



**Simplyread™**

**BRAND** 20 years of IVD experience; intellectual property rights;

### Technology and Innovation

- Reagent R&D Center in US;
- Instrument technology from Freiburg Germany;
- 52 international patents and 12 invention patents.

### Facility:

World-class high capacity manufacturing facility ( 10000 sq.m.) with 1500 sq.m. current GMP plants and strict quality control system.

### Production capacity:



- 150000 Tests/day;
- Certified by TUV Rheinland with ISO 13485 and CE marks.



### GLOBAL SERVICE NETWORK



### TEST LIST



R&D: 3938 Trust Way, Hayward, CA 94545, USA  
Production base: No.3 Building,Fangyi Tech-park, 313 zhuijiang Rd, Shijiazhuang city,P.R.China  
Service Hotline: 400-0191-606

Tel: +86 311 83832628 | Fax: +86 311 83859087  
Email: IMD@hipro.us | Website: www.hipro.us

# A1 [ INNOVATIVE DESIGN ]

## Automated Immunoassay System



**labONE™**

- Portable ·
- Multi methodologies ·
- Rapid quantitative results ·
- Local language selection ·
- Build-in Battery(optional) ·

**Hipro®**

Hipro Biotechnology Co.,Ltd



Originate in Silicon Valley