



PRODUCT NAME	TARGET COMPOUNDS	CUTOFF LEVELS(ng/ml)
THC (MARIJUANA)	11-NOR- $\Delta^9$ -THC-9-COOH	50
OPI (OPIATES)	MORPHINE	300/2000
COC (COCAINE)	BENZOYLECGONINE	300
AMP	D-AMPHETAMINE	1000
mAMP	D-METHAMOHETAMINE	1000
MDMA (ECSTASY)	METHYLENEDIOXYMETHAMPHETAMINE	500
MDN	METHADONE	300
BZO	OXAZEPAM	300
BAR (BARBITURATES)	SECOBARBITAL	300
TCA (TRI-CYCLIC ANTIDEPRESSANTS)	NORTRIPTYLINE	1000
PCP	PHENCYCLIDINE	25
PPX (PROPOXYPHENE)	NOR-PROXYPHENE	300
BUP (BUPRENORPHINE)	NORBUPRENOPHINE	10

## Nano-Check™ DAT Ordering Information

Drugs of Abuse Tests					
Catalog No.	Product Name	Specimen	Target Drug	Format	Test/Box
NDAD4183	Nano-Check™ Rapid OPI	Urine	Morphine (OPI)	Cassette/Strip	25/50
NDAD4283	Nano-Check™ Rapid COC	Urine	Cocaine (COC)	Cassette/Strip	25/50
NDAD4383	Nano-Check™ Rapid AMP	Urine	Amphetamine (AMP)	Cassette/Strip	25/50
NDAD4483	Nano-Check™ Rapid MET	Urine	Methamphetamine (MET)	Cassette/Strip	25/50
NDAD4583	Nano-Check™ Rapid THC	Urine	Marijuana (THC)	Cassette/Strip	25/50
NDAD4683	Nano-Check™ Rapid PCP	Urine	Phencyclidine (PCP)	Cassette/Strip	25/50
NDAD4783	Nano-Check™ Rapid TCA	Urine	Tricyclides (TCA)	Cassette/Strip	25/50
NDAD4883	Nano-Check™ Rapid BZO	Urine	Benzodiazepine (BZO)	Cassette/Strip	25/50
NDAD4983	Nano-Check™ Rapid BAR	Urine	Barbiturate (BAR)	Cassette/Strip	25/50
NDAD5083	Nano-Check™ Rapid MDN	Urine	Methadone (MDN)	Cassette/Strip	25/50
NDAD5183	Nano-Check™ Rapid EXT	Urine	Ecstasy (MDMA)	Cassette/Strip	25/50
NDAD5184	Nano-Check™ Rapid PPX	Urine	Propoxyphene (PPX)	Cassette/Strip	25/50
NDAD5185	Nano-Check™ Rapid BUP	Urine	Buprenorphine (BUP)	Cassette/Strip	25/50
NDAD200U	Nano-Check™ Rapid DAT2	Urine	Any 2 Drugs Combination Panel	Cassette	25
NDAD300U	Nano-Check™ Rapid DAT3	Urine	Any 3 Drugs Combination Panel	Cassette	25
NDAD400U	Nano-Check™ Rapid DAT4	Urine	Any 4 Drugs Combination Panel	Cassette	25
NDAD501U	Nano-Check™ Rapid DAT5M	Urine	THC/COC/OPI/MET/PCP	Cassette	25
NDAD504U	Nano-Check™ Rapid DAT5B	Urine	BAR/BZO/MDN/AMP/TCA	Cassette	25
NDAD500U	Nano-Check™ Rapid DAT5U	Urine	Any 5 Drugs Combination Panel	Cassette	25
*NDAD601U	Nano-Check™ Rapid DAT6M	Urine	THC/COC/OPI/BAR/BZO/AMP	Cassette	25
*NDAD801U	Nano-Check™ Rapid DAT8U	Urine	Any 8 Drugs Combination Panel	Cassette	25
*NDAD101U	Nano-Check™ Rapid DAT10U	Urine	Any 10 Drugs Combination Panel	Cassette	25

- Any combination of two to five drugs of abuse listed in this table is available on the strip
- Nano-Ditech can provide multi or single line drug tests in three different formats : multi cassettes, dipstick strips, single cassette, and also in your own device.