

# Nano-Check™

# DAT Drugs of Abuse Panels

## Applications

The Nano-Check™ DOA test is a One-step immunochromatography assay that detects the presence of any of 13 different illicit drugs in urine samples within 10 minutes. The Nano-Check™ DOA test is the most comprehensive multidrug panel format in the DOA test industry. The product is easy to use and can screen as many as 10 different drugs of abuse simultaneously.

## Accuracy and Quality Consistency

- Rigorous quality Control System
- Extensively checked cross reactivity to more than 500 chemical substances
- Clear cut discrimination above and below cutoff value
- Manufactured by Nano-Ditech Corp. in USA

## Convenient

- Simple One step Assay
- Small application volume of samples
- One strip for up to five classes of abused drugs
- Read by portable reader or eye vision
- Easy to use Nano-Checker 710N Reader, results are automatically determined and displayed or printed as 'Positive' or 'Negative'
- Clear visual signal, minimizes reading error

## Economic

- No calibration, built in control
- No maintenance cost
- Multi test in one cassette
- Variable formats for up to 11 tests

