



Third Party
Quality Controls



External Quality
Assessment



Interlaboratory
Data Management



Molecular
IQC & EQA



ACUSERA QUALITY CONTROLS



Immunoassay Premium Plus Quality Controls

The **Randox Acusera Immunoassay Premium Plus** quality control covers an impressive 55 analytes with assayed instrument specific target values and ranges provided for up to 49 analytes.

The unique combination of routine tumour markers, therapeutic drugs and Vitamin D allows laboratories to dramatically reduce the number of controls required while ultimately keeping costs to a minimum.

Three levels of control are available with analytes present at clinically significant levels.

Description	Cat. No.	Size	Analytes
Immunoassay Premium Plus Level 1	IA3109	12 x 5 ml	54
Immunoassay Premium Plus Level 2	IA3110	12 x 5 ml	54
Immunoassay Premium Plus Level 3	IA3111	12 x 5 ml	54

[Learn More](#)

FEATURES & BENEFITS



Lyophilised control format ensuring enhanced stability.



Manufactured using 100% human serum, this control mimics the performance of the patient sample.



Assayed target values are provided including Ferritin and Vitamin B12 levels suitable for anaemia monitoring and an ultra-low TSH in the L1 Control.



Control contains routinely run tumour markers including AFP, CA15-3, CA19-9, CA-125, PSA and Free PSA, providing effective consolidation.



Ferritin and Vitamin B12 levels suitable for anaemia monitoring

[Learn More](#)

[Contact Us](#)

RESOURCE HUB

The full extensive range of our QC literature material in the one place.



Randox Laboratories Ltd.

55 Diamond Road,
Crumlin,
County Antrim,
BT29 4QY,
United Kingdom

Tel: +44 (0) 28 9442 2413

Website: www.randox.com

Email: market@randox.com

Copyright © Randox Laboratories Ltd, All rights reserved.

Want to change how you receive these emails? [Update your preferences](#) or [unsubscribe from this list](#)