



Third Party  
Quality Controls



External Quality  
Assessment



Interlaboratory  
Data Management



Molecular  
IQC & EQA



## SERUM INDICES CONTROL



The **Randox Acusera Serum Indices control** is designed to be used to monitor an IVD instrument's response in the detection of HIL samples. This control can be utilized in laboratory interference testing to assist in improving error detection of pre-analytical errors affecting clinical chemistry testing.

Description	Pack Size	Cat. Number
Serum Indices Control	4 x 5 ml	SI10448

## FEATURES & BENEFITS

- **All HIL samples are provided in one control pack**, helping labs to consolidate testing and remove the need to purchase multiple kits
- **Lyophilised** for enhanced stability
- Reconstituted stability of **14 days at 2°C – 8°C** helping to minimise waste
- Superior manufacturing using **human serum**, ensures a commutable quality control is delivered
- **2-year shelf life** from date of manufacture



[Learn More](#)

[Contact Us](#)



## RELATED PRODUCTS

### RIQAS



The **RIQAS Serum Indices EQA programme** is designed for the pre-analytical assessment of Haemolytic, Icteric and Lipemic (HIL) interferences. It gives laboratories the option to report either quantitative or semi-quantitative results for HIL parameters whilst providing an assessment on how these interferences can impact on up to 25 routine chemistry parameters.

ALP | ALT | AST | D. Bilirubin | T. Bilirubin | Calcium | Chloride | Cholesterol | CK NAC  
| Creatinine | GGT | Glucose | HDL | Iron | Lactate | LDH | Lipase | Magnesium |  
Phosphate | Potassium | Protein Total | Sodium | Triglycerides | Urea | Uric Acid

Description	Frequency	Size	Cat. No
Serum Indices	Bi-monthly <i>(3 samples per challenge)</i>	2 x (9x1ml)	RQ9194

## FEATURES & BENEFITS

- **Lyophilised** for enhanced stability
- **Human based serum** ensuring commutable sample matrix
- **Bi-monthly reporting**
- HIL parameters include the option of **quantitative or semi-quantitative** reporting
- **Submit results and view reports online** via RIQAS.Net
- **Rapid report turnaround** ensures corrective actions can be taken swiftly

[Contact Us](#)

[Download Flyer](#)

**RANDOX**

**RIQAS** 

**RIQAS ROADSHOWS  
UK & IRELAND  
SEPTEMBER 2024**



Presented by:  
**Stephen Doherty**  
RIQAS Technical Expert



**Randox Laboratories Ltd.**

55 Diamond Road,  
Crumlin,  
County Antrim,  
BT29 4QY,  
United Kingdom

**Tel:** +44 (0) 28 9442 2413

**Email:** [qualitycontrol@randox.com](mailto:qualitycontrol@randox.com)

**Web:** <https://www.randox.com/>

Want to change how you receive these emails?

You can [update your preferences](#) or [unsubscribe](#)